

**In the Specification:**

Please replace the paragraph beginning at page 7, line 14, with the following:

--Figure 2 provides a schematic representation of the differentially methylated region (DMR) of the human *IGF2-H19* region. The two 450-bp repeat (A1 and A2) and seven 400-bp repeat (B1-B7) units are shown. The potential methylation sites on the upper strand DNA of the studied region are represented by open circles. The studied single nucleotide polymorphism (SNP) site (A/G) is indicated by an open box. Open arrows represent the location of the forward (for) and reverse (rev) primers in PCR reactions specific for the methylated (M) and unmethylated (U) alleles, respectively. Sequences of these MSP primers are shown (SEQ ID NOS:7-10). Sequence differences between bisulfite-treated DNA and untreated DNA are highlighted in bold italics and sequence differences between methylated (paternally-inherited) and unmethylated (maternally-inherited) DNA are underlined in bold.--

Please replace the paragraph beginning at page 15, line 22, with the following:

--MSP assays were modified from the protocol as described by Herman *et al, supra*. The primers M-for (5'-GCGAGCGTAGTATTTTCGGC-3'; SEQ ID NO:1) and M-rev (5'-AACCAAATAACCTATAAACCTCTACG-3'; SEQ ID NO:2) were designed for the methylated sequence, while the primers U-for (5'-GTTGTGAGTGTAGTATTTTGTT-3'; SEQ ID NO:3) and U-rev (5'-CAAATAACCTATAAACCTCTACA-3'; SEQ ID NO:4) were designed for the unmethylated sequence. Five  $\mu$ l bisulfite-treated DNA was added to a 50  $\mu$ l PCR reaction containing 5  $\mu$ l 10x TaqMan buffer A (PE Applied Biosystems), 2 mM MgCl<sub>2</sub>, 10 pmol dNTPs, 20 pmol each of the corresponding MSP primers and 1.25 U AmpliTaq Gold DNA polymerase (PE Applied Biosystems). Reaction mixtures were thermal cycled

(methylated allele: 95 °C for 45 sec, 58 °C for 30 sec, 72 °C for 20 sec; unmethylated allele: 95 °C for 45 sec, 50 °C for 30 sec, 72 °C for 20 sec) for 45 cycles, with an initial denaturing step of 8 min at 95 °C. PCR products were then analyzed by agarose gel electrophoresis.--

Please replace the paragraph beginning at page 17, line 32, with the following:

--The DMR in the human IGF2-H19 locus contains two 450-bp repeat and seven 400-bp repeat units (Nakagawa et al., *Proc Natl Acad Sci U S A* **98**:591-596 (2001)) (Fig. 2). An A/G SNP within the DMR (Nakagawa et al., *supra*) was selected as a marker in our investigation (Fig. 2). Polymerase chain reaction (PCR) was used to amplify the SNP in both maternal and fetal DNA samples. Primers were designed using the sequence of the *Homo sapiens H19* gene (Genbank accession number AF125183). Typically, 2 to 5 µl eluted DNA, purified from maternal buffy coat, cord buffy coat or amniotic fluid was added to a 25 µl PCR reaction containing 2.5 µl 10x TaqMan buffer A (PE Applied Biosystems), 3 mM MgCl<sub>2</sub>, 6.26 pmol dNTPs, 5 pmol primers (forward: 5'-ggACGAATTGGTTAGTT-3' (SEQ ID NO:5); reverse: 5'-AGGCAATTGTCAGTTCAGTAA-3' (SEQ ID NO:6)) and 0.625 U AmpliTaq Gold DNA polymerase (PE Applied Biosystems) (95 °C for 8 min followed by 35 cycles of 95 °C for 1 min, 56 °C for 20 sec, 72 °C for 20 sec). For the forward primer, the nucleotides in upper case corresponded to positions 7927 to 7944 of the *H19* sequence (Genbank accession number AF125183). For the reverse primer, the nucleotides were complementary to positions 8309 to 8329 of the *H19* sequence. PCR products were then analysed by agarose gel electrophoresis and DNA sequencing.--

Please replace the paragraph beginning at page 18, line 26, with the following:

--MSP assays were modified from the protocol as described (Herman et al. 1996). Five µl bisulfite-treated DNA was added to a 50 µl PCR reaction containing 5 µl

10X TaqMan buffer A (PE Applied Biosystems), 2.5 mM MgCl<sub>2</sub>, 10 pmol dNTPs, 20 pmol each of the corresponding MSP primers (Fig. 2) and 1.25 U AmpliTaq Gold DNA polymerase (PE Applied Biosystems). The primers M-for and M-rev (Fig. 2; SEQ ID NOS:7 and 8) were designed for the methylated sequence, while the primers U-for and U-rev (Fig. 2; SEQ ID NOS:9 and 10) were designed for the unmethylated sequence. Reaction mixtures were thermal cycled (methylated allele: 95 °C for 45 sec, 55 °C for 20 sec, 72 °C for 20 sec; unmethylated allele: 95°C for 45 sec, 49 °C for 20 sec, 72 °C for 20 sec) for 50 (buffy coat and amniotic fluid DNA) or 56 (plasma DNA) cycles, with an initial denaturing step of 8 min at 95 °C. PCR products were then analyzed by agarose gel electrophoresis. Reaction products were purified using Microspin S-300 HR columns (Amersham Pharmacia) for DNA sequencing or the primer extension assay.--

Please replace the paragraph beginning at page 19, line 11, with the following:

--Two  $\mu$ l of the purified MSP product was added to a 25  $\mu$ l reaction containing 50  $\mu$ M ddATP (2',3'-dideoxyadenine triphosphate), 50  $\mu$ M dGTP, 50  $\mu$ M dTTP, 0.2 pmol Cys-5-labeled primer (5'-GGGTTATTGGGAATAGGATATT-3'; SEQ ID NO:11), 4 U Thermo Sequenase (Amersham Pharmacia) and 1.43  $\mu$ l concentrated buffer. Reactions were thermal cycled for 40 cycles (95°C for 30 sec, 51°C for 20 sec, 72°C for 20 sec). The Cys-5-labeled primer was 25 nucleotides (nt) in length and the polymorphic site was 2 nt away from the 3'-end of the primer. For the A allele, the incorporation of the ddATP at this polymorphic site would produce chain termination, thus resulting in an extension product of 27 nt (i.e., 25+2 nt). For the G allele, chain extension would continue until the next A residue which was 5 nt away from the 3'-end of the primer, thus resulting in an extension product of 30 nt (i.e., 25+5 nt). Reaction products were electrophoresed using a 14% denaturing polyacrylamide gel and analysed using an ALF Express Sequencer (Amersham Pharmacia). Data were analysed by the AlleleLinks program (Amersham Pharmacia).--